



REFERENCES

Created: 8/17/20
Department of Emergency Medicine

Cardiac stabilization
3g calcium gluconate IV or 1g calcium chloride IV
- Can repeat q5-10 min x 3 total doses
- Obtain repeat Ca levels

Elimination
Kaliuresis - see algorithm above
Consider: sodium zirconium cyclosilicate 10 g PO
patiomer 8.4 g PO
or kayexlate 15 g PO
Dialysis - consult renal fellow

Intracellular shift
10 units insulin regular IV and 20mg albuterol neb
* if glucose < 250, give 2 amps of D50
(May prefer this choice if rapid effect is desired)

20 units insulin regular IV infusion + D10 infusion and
20mg albuterol neb
(May prefer this option if more gradual effect is desired)
POC glucose q30-60 min